

Product Features

Product Name : Pure Pro

Product Category : Skin Booster

Ingredients:

Polynucleotide - 22mg/ml

Hyaluronic Acid - 5.5mg/ml

Peptides

- Palmitoyl Pentapeptide-4

- Octapeptide-4

- Acetyl Hexapeptide-8

Amino Acids

- Glycine, Proline, L-Hydroxyproline

- Leucine, Glutamin, L-Lysine, Methionine

Package:

2.2ml x 1 syringe / box

Manufactured in South Korea



Benefits of PN

SKIN REPAIR

Boosts cell renewal and stimulates new cell growth, promoting collagen production.

HYDRATION

Works by attracting and retaining moisture in the skin, improving hydration.

REJUVENATION

Pure Pro

Polynucleotide + Hyaluronic Acid + Peptides

High grade PN Booster

Promotes a more youthful and rejuvenated complexion by boosting collagen levels.

ELASTICITY I M PROVE M E NT


Stimulates elastin to increase skin elasticity.

SKIN TONE I M PROVE M E NT

Improves pigmentation and skin tone by increasing the thickness and density of the dermis.

PORE TIGHTENING

Provides pore tightening effects by filling open pores from within the skin.






Hyaluronic Acid:

When hyaluronic acid is injected directly into the dermis, it immediately enhances the skin's moisturizing power, improving weakened elasticity and rough texture.

As a result, the skin becomes smoother and more radiant.

It can be safely combined with PN-based skin boosters, which are effective for regeneration, leading to optimal treatment results.





Peptides:

- •Acetyl Hexapeptide-8: A neurotransmitter-inhibiting peptide that effectively reduces wrinkles by acting on the skin's surface layer.
- •Octapeptide-4: Protects the skin from free radicals that cause aging, helping prevent oxidative damage and promoting a more vibrant complexion.
- •Palmitoyl Pentapeptide-4: Stimulates collagen production, elastic fibers, and hyaluronic acid, improving the skin's moisture content and retention. It also increases skin thickness and reduces fine lines.

Amino Acids:

PURE PRO contains key amino acids that enhance collagen synthesis.

INGREDIENTS	EFFECTS
Glycine Proline L-Hydroxyproline	Collagen production Growth factor generation Promotion of collagen synthesis
Leucine Glutamin L-Lysine Methionine	Skin fibroblast activation + Natural moisturizing factor





Procedure for Pure Pro

Method: Micro-Droplet Technique

Needle Size: 32G 4mm

Cannula Size: 22/23g - 50-70mm

Injection Areas:

Forehead

· Crow's feet

Cheek

Neck

· Back of hand

Elbow

How to Use:

• At 1 cm intervals, inject 0.025 ml to 0.05 ml per point.

• Inject a total of 2 ml throughout the entire procedure area.

Recommended Sessions 3 sessions (2 - 3 weeks interval)

How Long Does It Last: 6 months

Test Results – Before and After

Quick Effect

3~5 days after the first procedure, the balance of oil and water in the skin is controlled and the skin texture on the face becomes soft and clear.

Long - Lasting Effect

The effect created by Pure Pro lasts more than 6 month, and over time, the ingredients stimulate collagen regeneration, which increases skin elasticity.

Natural Result

It promotes collagen production in the body, leading to gradual and natural changes without filling in artificial substances.







